



**TECHNISCHE
UNIVERSITÄT
DRESDEN**

Faculty of Electrical and Computer Engineering, Communications Laboratory

Chair for RF and Photonics Engineering

Electrode-based implanted HBC channel characterization and SAR analysis

Q. Wang, X. Du, T. Bauer, M. Baerhold, D. Plettemeier

Chair for RF and Photonics Engineering, Communication Laboratory

TU Dresden

URSI GASS 2020




**DRESDEN
concept**
Exzellenz aus
Wissenschaft
und Kultur

HBC (Human Body Communications):


- a novel rising transmission technique
- using the conductive properties of the body for electrical signal transfer
- small power requirement and absence of antenna
- potential alternative scheme of wireless RF communication

Electrode-based implanted HBC:

- pair of electrodes used
- electric field propagating through the body via the galvanic coupling
- to explore the transmission channel gain and the SAR analysis
- frequency range of interest of 5MHz-50MHz
- quasi-electrostatic approximation




Concept of transmission using the human body as electrical transmission channel




Electrode pill models


- m1: $D=0.6\text{cm}$, $L_1=5\text{cm}$ (big pill), $L_2=0.8\text{cm}$ (long electrode)
- m2: $D=0.6\text{cm}$, $L_1=5\text{cm}$ (big pill), $L_2=0.2\text{cm}$ (short electrode)
- m3: $D=0.6\text{cm}$, $L_1=2\text{cm}$ (small pill), $L_2=0.8\text{cm}$ (long electrode)
- m4: $D=0.6\text{cm}$, $L_1=2\text{cm}$ (small pill), $L_2=0.2\text{cm}$ (short electrode)




m1



m2




m3



m4

Four fabricated electrode prototypes



Electric potential distribution for the implanted electrode with low frequency time domain Electro-Quasi-Static solver in CST



Phantom measurement with balanced Balun PI100



Water	52.4%
NaCl	1.4%
Sugar	45%
Sodium benzoate	1.1%

phantom recipe in HBC band

$$Gain[dB] = 20 * \log_{10} \frac{|U_{RX1} - U_{RX2}|}{|U_{TX1} - U_{TX2}|}$$


Channel gain characterization

Measured channel gain in phantom for electrode pair m1



Relative good matching impedance characteristics inside the phantom
S11 are below -10dB


- Frequency independent gain characteristics.
- Same conclusion drawn as well for m2, m3 and m4.
- Decreasing gain with increasing transmission distance.



m1 (long rod and long electrode)
m2 (long rod and short electrode)

- Increasing the electrode length by 0.6cm (electrode length difference of m1 and m2) brings a channel gain enhancement of approximately 10dB.



Inter-electrode distance influence on the channel gain




m1 (long rod and long electrode)
m3 (short rod and long electrode)

- Increasing the inter-electrode distance by 3cm (difference of m1 and m3) brings a channel gain enhancement of approximately 10dB.

Influence of the electrode diameter variation




- Increasing the electrode diameter by 0.6cm brings a channel gain enhancement of approximately 8dB.



m1 (long rod and long electrode)
m4 (short rod and short electrode)


- Measured channel gain shows frequency independence (10MHz-50MHz).
- At 10cm, 15dB PL difference can be observed of big and small electrodes.
- For m1 electrode, at 10cm, the PL is around 45dB.
- For m4 electrode, at 10cm, the PL is around 60dB.

Measured channel gains in phantom with electrode pair distance 5cm




- The larger electrode m1 demonstrate the relative best channel gain characteristics.

Comparison of alignment and 45 degree non-alignment, with 7 cm electrode pair distance




- Compared with alignment, 45 degree displacement corresponds to approximately 1.5dB@m1, 2.5dB@m2, 2.5dB@m3 additional loss.
- The smaller the pill, the larger the additional loss variation. (m3 and m4)
- The larger the metal electrode, the smaller the additional loss. (m1 vs. m2)
- Electrode vertical displacement lose the transmission connection.

m1 (long rod and long electrode)




Homogeneous CAD torso model
(with 2/3 muscle dielectric properties)



10-g SAR distribution on cross-section
cut plane of the voxel torso model


- 0.2V potential excitation at 21MHz
- Electrode pairs locating in muscle tissue region

10-g averaged SAR peak comparison of homogeneous CAD model and heterogeneous Voxel model (with 0.2V potential excitation)



- SAR with homogeneous model (2/3 muscle dielectric properties) and heterogeneous model show similar results.
- SAR in mW/kg, much smaller than the safety guideline 2 W/kg.

10-g averaged SAR peak comparison of four electrode prototypes with 1 mA current excitation



- Smaller metal electrode (m2 and m4) demonstrates higher SAR, since the resulting E-field is higher around the electrode.
- SAR highly localized, rod length make less influence on the SAR (m1 vs. m3 and m2 vs. m4).
- Larger implanted electrode demonstrates the relative better channel gain as well as lower localized SAR peak.